

Vbeam® Prima

Technology You Trust.
Innovation You Need.




 **CANDELA™**
Science. Results. Trust.




Trustworthy technology.
Thoughtfully designed for
vascular treatment and more.

Versatile treatments

The Vbeam® Prima device is a highly efficient pulsed-dye laser (PDL) used by healthcare providers all over the world to treat a variety of indications for both face and body, including benign vascular and pigmented lesions, with high patient tolerability and a low incidence of side effects.¹⁻¹²



A legacy of innovation



Two wavelengths, more
treatment parameters

Proprietary 595 nm wavelength

The 595 nm wavelength deeply penetrates the skin to reach targeted blood vessels.⁴ Its energy is absorbed by oxyhemoglobin to coagulate and clear vessels with more tolerability and fewer instances of melanin absorption.¹³



8 micropulses

8 micropulses allow the physician to treat with or without purpura, increasing treatment success or reducing patient downtime when required.

Additional 1064 nm wavelength


Treat deep blue leg veins, venous lakes, and wrinkles¹²

- Treats blue veins across face and body
- Minimises the appearance of wrinkles



New features that benefit both the provider and patient experience

Achieve greater results in less time



Larger max treatment spot size and greater max energy with 595 nm

- Faster, more efficient treatments
- A treatment that requires 100 pulses with a 12 mm spot size can require as few as 64 pulses with a 15 mm spot size


	VBEAM® PERFECTA	VBEAM® PRIMA
Maximum energy	8 J	12 J
Maximum spot size	12 mm	15 mm
Area covered	1.13 cm²	1.77 cm²
Fluence at 12mm	7.00 J/cm²	9.75 J/cm²
Fluence at 15mm	-	6.75 J/cm²



EverCOOL

Cryogen-based Dynamic Cooling Device™ (DCD™) and EverCool™ contact cooling options

- Two types of cooling for maximum versatility and epidermal protection
- DCD™ scales with fluence to automatically administer consistent epidermal protection



Smart dye life management

- Dye life meter provides data on remaining dye levels
- Avoid unexpected treatment interruptions from dye loss
- Allows user to schedule preventative maintenance



Once-daily system calibration per wavelength^b

Device is patient and treatment ready for the entire day, reducing downtime between patients; no need to recalibrate when changing spot sizes or fluences



Zoom handpiece

A versatile handheld applicator that allows for very targeted spot size adjustments in increments as small as 0.5 mm



Wi-Fi connectivity

Accurate remote service diagnostics for faster equipment servicing



Guided user interface

Intuitive, easy-to-use software with quick access to saved, favourite treatment settings




^bAuto-calibration on start up.


Powerful results begin with Vbeam® PrimaPlatform

The Vbeam® Prima Platform is the trusted, proven PDL across indications¹⁴

Telangiectasia




Before
3 months post Vbeam Prima treatment.
Photos courtesy of E. Victor Ross, MD.



After*
3 months post Vbeam Prima treatment.
Photos courtesy of E. Victor Ross, MD.

Pigmentation



Before
3 months post Vbeam Prima treatment.
Photos courtesy of E. Victor Ross, MD.



After*
3 months post Vbeam Prima treatment.
Photos courtesy of E. Victor Ross, MD.



Diffuse Redness



Before
3 months post Vbeam Prima treatment.
Photos courtesy of E. Victor Ross, MD.




After*
3 months post Vbeam Prima treatment.
Photos courtesy of E. Victor Ross, MD.

*All photos are unretouched. Individual results may vary.

Demand the best for your patients

An estimated 16 million people in the United States suffer from rosacea¹⁵



Significant improvements in rosacea symptoms with Vbeam® Prima Platform

In a clinical study, rosacea symptoms were significantly reduced from moderate to mild.^{17,18}

- ✓ Difficult-to-treat erythema was reduced from severe to mild
- ✓ 84% of patients had >40% rosacea improvement



Minimal adverse side effects reported¹⁷

On a scale of 1 to 10 (10=max), patients reported an average pain score of 5.6 ± 1.8 . Most patients reported side effects such as mild edema, mild to moderate erythema, and mild to moderate purpura. These effects resolved within a few days to a week after treatment without intervention.



Before



After*

2 months post four treatments of Vbeam Prima.
Photos courtesy of Eric F. Bernstein, MD, MSE.

*All photos are unretouched. Individual results may vary.



SYSTEM SPECIFICATIONS ^c	
Laser wavelengths	PDL: 595 nm Nd:YAG laser: 1064 nm
Laser pulse repetition rate	PDL: Up to 1.5 Hz Nd:YAG laser: Up to 10 Hz
Laser pulse duration	PDL: 0.45-40 ms 8 micropulses technology Nd:YAG laser: 0.5-60 ms
Maximum energy	PDL: 12 J Nd:YAG laser: 45 J
Method of optical output	Lens-coupled optical fiber with user-selectable spot sizes
Networking method	Wi-Fi
Dimensions (H x W x D)	53 x 20 x 33 in / 135 x 51 x 84 cm
Weight	280 lbs / 127 kg
Electrical requirements	220-240 VAC, 24 A max, 50/60 Hz, single phase
Patented Dynamic Cooling Device (DCD™) Integrated controls, cryogen container, and handpiece with distance gauge	
Cryogen	HFC 134a
DCD spray duration	User-adjustable range: Off, 20-100 ms
DCD delay duration	User-adjustable range: 10-150 ms
DCD post-spray duration	User-adjustable range: Off, 10-50 ms
Beam spot sizes	1-5, 3-15 mm Zoom handpiece 3 x 10 mm elliptical spot attachment
EverCool Contact Cooling User-controlled, adjustable sapphire cooling tip for use pre, during, and post pulse	
Temperature range	10°C-20°C
Beam spot sizes	3-15 mm Zoom handpiece

- Large 15 mm spot size and up to 50% greater power at 12 J with 595 nm - for faster, more efficient treatments
- Once-daily system calibration per wavelength^b - for less downtime between treatments, and more efficient treatments when using a variety of settings
- Smart dye life management - for real-time data on remaining dye levels
- 1064 nm wavelength - for deep blue vessel clearance and wrinkle reduction
- Cryogen-based DCD™ and EverCool contact cooling options - for maximum epidermal protection and treatment versatility

For more information about how the Vbeam Prima system may help achieve your practice goals, contact your local Candela sales professional or visit candelamedical.com

^bAuto-calibration on start up.

^cSystem specifications subject to change without notice.

References. **1.** Vbeam 510(k) clearance (K183452), January 2004. **2.** Vbeam 510(k) clearance for pigmented lesion handpiece accessory (K183452), July 2005. **3.** Bernstein EF, Kligman A. Rosacea treatment using the new-generation, high-energy, 595 nm, long pulse-duration pulsed-dye laser. *Lasers Surg Med.* 2008;40(4):233-239. **4.** Woo SH, Ahn HH, Kim SN, Kye YC. Treatment of vascular skin lesions with the variable-pulse 595 nm pulsed dye laser. *Dermatol Surg.* 2006;32(1):41-48. **5.** Chapas AM, Eickhorst K, Geronemus RG. Efficacy of early treatment of facial port wine stains in newborns: a review of 49 cases. *Lasers Surg Med.* 2007;39(7):563-568. **6.** Jasim ZF, Woo WK, Handley JM. Long-pulsed (6-ms) pulsed dye laser treatment of rosacea-associated telangiectasia using subpurpuric clinical threshold. *Dermatol Surg.* 2004;30(1):37-40. **7.** Jørgensen GF, Hedelund L, Hædersdal M. Long-pulsed dye laser versus intense pulsed light for photodamaged skin: a randomized split-face trial with blinded response evaluation. *Lasers Surg Med.* 2008;40(5):293-299. **8.** Halachmi S, Israeli H, Ben-Amitai D, Lapidot M. Treatment of the skin manifestations of hereditary hemorrhagic telangiectasia with pulsed dye laser. *Lasers Med Sci.* 2014;29(1):321-324. **9.** Yu W, Ma G, Qiu Y, et al. Prospective comparison treatment of 595-nm pulsed-dye lasers for virgin port-wine stain. *Br J Dermatol.* 2015;172(3):684-691. **10.** Galeckas KJ, Ross EV, Uebelhoefer NS. A pulsed dye laser with a 10-mm beam diameter and a pigmented lesion window for purpura-free photorejuvenation. *Dermatol Surg.* 2007;34:1-6. **11.** Madan V, Ferguson J. Using the ultra-long pulse width pulsed dye laser and elliptical spot to treat resistant nasal telangiectasia. *Lasers Med Sci.* 2010;25(1):151-154. **12.** Indications for 1064 nm wavelength. Candela, data on file. **13.** Bernstein EF. The pulsed-dye laser for treatment of cutaneous conditions. *G Ital Dermatol Venereol.* 2009;144(5):557-572. **14.** Ross EV. Vbeam Prima before and after photos. Candela, data on file. **15.** National Rosacea Society website. <https://www.rosacea.org/patients/index.php>. Accessed November 8, 2017. **16.** Consumer rosacea laser attitude and behavior exploratory: final report. September 7, 2017. BuzzBack Market Research. **17.** Bernstein EF, Schomacker K, Paranjape A, Jones CJ. Pulsed dye laser treatment of rosacea using a novel 15 mm diameter treatment beam. *Lasers Surg Med.* 2018;doi:10.1002/lsm.22819. **18.** Bernstein EF, Schomacker K, Paranjape A, Jones CJ. Pulsed dye laser treatment of rosacea using a novel 15 mm spot size. Candela, data on file.