

PICK YOUR POWERHOUSE

DURABILITY. VERSATILITY. RESULTS.



The Gentle Pro Series™

PICK YOUR POWERHOUSE

DURABILITY. VERSATILITY. RESULTS.



For over **25 years**, aesthetic practices have successfully built their business with help from the power and performance of the Gentle series of lasers.¹⁻³

Whether you're seeking an entry level hair removal laser that allows you to upgrade as your business grows, or a best-in-class all-in-one device from the start, the Gentle Pro Series™ has a powerhouse device that's right for you.

The Industry Choice

The Gentle series of devices are well established as the aesthetic industry's best-in-class laser hair removal systems:

- **The world's leading brand of hair removal lasers⁴:** Time-tested for over 25 years³
- **Proven results:** Significant hair reduction after just 3 treatments¹⁻³

The Gentle Pro Series™

PICK YOUR POWERHOUSE

The Gentle Pro Series™

GentleMax Pro Plus® | GentleMax Pro® | GentleLase Pro-U™ | GentleYag Pro-U™



+ More results for your patients

- **More power:** Enables higher fluence to efficiently remove hair, from the coarsest to the finest, in less time.¹⁻³
- **More speed:** Treats more patients in less time with larger spot sizes and higher repetition rates for quicker treatments.¹⁻³
- **More spot sizes:** Spot sizes range from 1.5 mm to 26 mm for customized treatments across conditions and body areas.¹⁻³
- **More fine hair removal:** Short pulse durations to improve targeting of fine, resistant hair.¹⁻³
- **More satisfaction, comfort, and long-lasting results:** Precision DCD™ and air cooling options provide comfortable treatments for high patient satisfaction.³

+ More possibilities for your practice

- While known for expert-level hair removal, together the 755 nm and 1064 nm wavelength lasers expand your practice offerings with treatment of: benign pigmented lesions, vascular lesions, redness and facial vessels, spider and leg veins, wrinkles, and onychomycosis (nail fungus).¹

The Gentle Pro Series™

PICK YOUR POWERHOUSE

MEET YOUR POWERHOUSE



GentleMax Pro Plus®

Candela's latest evolution in the series is built on an exceptional multi-technology platform that provides a best-in-class solution, with treatment versatility for in-demand procedures. Attract and retain patients, all while adding additional revenue streams with a high ROI.³

The GentleMax Pro Plus device includes all of the GentleMax Pro® features, plus:

- More energy in larger spot sizes.¹⁻³
- Unprecedented speed with proven results.¹⁻³
- Shorter pulse durations to treat fine, residual hair.¹⁻³
- Faster repetition rates for improved speed.¹⁻³



GentleMax Pro®

Combining the power and treatment capabilities of the GentleLASE Pro-U™ and GentleYAG Pro-U™ devices, this trusted dynamo includes both 755 nm and 1064 nm wavelengths for treatment versatility.¹⁻³

- Offers a full range of treatment capabilities beyond hair removal, including benign pigmented lesions, vascular lesions, wrinkles, telangiectasia, and onychomycosis.¹
- Treats all skin types.¹

The Gentle Pro Series™

PICK YOUR POWERHOUSE



GentleLase Pro-U™

This stand-alone 755 nm powerhouse is an ideal solution for those starting a practice or new to laser hair removal.

- Extends treatments beyond hair removal to include benign pigmented or vascular lesions.
- Optimal for lighter skin tones.¹
- Upgradeable for superior performance and additional clinical indications^{2,3} - add a 1064 nm laser when ready.



GentleYag Pro-U™

A stand-alone Nd:YAG 1064 nm laser platform that achieves successful treatment of the most sought-after cosmetic conditions.

- An ideal solution for hair removal, vascular lesions, wrinkles and diffuse redness.¹⁻³
- Optimal for darker skin tones.¹
- Upgradeable – add a 755 nm wavelength as your practice grows.²⁻³



The Gentle Pro Series™

PICK YOUR POWERHOUSE

FOUR POWERHOUSE LASERS FOR PRACTICE GROWTH²⁻³

THE GOLD STANDARD FOR HAIR REMOVAL.
PLUS EXPANDED TREATMENT CAPABILITIES.

	GentleMax Pro Plus [®]	GentleMax Pro [®]	GentleLase Pro-U [™]	GentleYAG Pro-U [™]
Ideal Practice Types	Busy Laser Hair Removal	Busy to Occasional Laser Hair Removal	Occasional Laser Hair Removal or Practices in Need of a Supplemental Device	
Wavelengths	755 nm, 1064 nm	755 nm, 1064 nm	755 nm	1064 nm
Treatment Versatility & Speed	++++	+++	++ Expandable	++ Expandable
Large range of spot sizes	1.5–26 mm	1.5–24 mm	3–24 mm	1.5–24 mm
Short pulse durations	0.25–100 ms	0.25–100 ms	3–100 ms	2.5–100 ms
External calibration	✓	✓	✓	✓
DCD [™] or ACC Cooling Options	✓	✓	✓	✓
Upgradable	--	--	✓	✓
Skin Types	All Skin Types	All Skin Types	All Skin Types Optimal for Fitzpatrick Skin Types I-III	All Skin Types Optimal for Fitzpatrick Skin Types IV-VI

The Gentle Pro Series[™]

GentleMax Pro Plus^{*} | GentleMax Pro^{*} | GentleLase Pro-U^{*} | GentleYag Pro-U^{*}

For more information about how the Gentle Series[™] devices can benefit your practice, contact your local Candela sales professional, or visit candelamedical.com.

1. Gentle Family , CE mark. 2. Gentle Family User Manuals: GentleMax Pro, GentleLase Pro-U, GentleYag Pro-U, LBL102517, GentleMax Pro Plus , LBL102519. 3. Candela Data on File, 2024. 4. Wizer Survey of Dermatologists.

