

eTwo System – A smart choice

- ① Intelligent Feedback System provides real time information on skin condition, impedance and the energy delivery of each pulse for consistency and reproducibility
- ② Minimal downtime so patients can quickly resume their normal activities
- ③ Safe and effective¹
- ④ Quick treatments: typically up to 1 hour depending on area treated in combined Sublative and Sublime protocol
- ⑤ Short treatment protocol: Visible improvement after the first session
- ⑥ Lightweight, compact device that easily transports between locations



eTwo® System Specification

Sublime Energy	Up to 200 J/cm ³
Sublative Energy	Up to 100 mJ/pin
Sublime Spot Size	12 x 8 mm
Sublative iD Spot Size	Standard 64 - 12 x 12 mm Focal 44 - 11 x 3 mm
Disposable Sublative iD	Standard 64 up to 400 pulses; Focal 44 up to 200 pulses
User Interface	15" touch screen
System Dimensions Table Top	53 x 30 x 40 cm (21.8" x 11.8" x 15.7")
System Weight	8.5 kg / 18.7 lbs
Electrical Requirements	100-240 VAC; 3A; 50-60 Hz

For more information about how the eTwo system may help achieve your practice goals, contact your local Candela sales professional or visit candelamedical.com

Facial rejuvenation in a compact solution



Developed with your patients in mind

The promise of noticeable results with little to no downtime and lower procedure cost is spurring interest in fractional methods of skin rejuvenation. Patients are increasingly requesting these treatments to improve their appearance. The eTwo system not only enables your practice to address this growing demand, but also adds non-invasive wrinkle treatment¹ to your revenue stream.

Meet the Growing Demand

The eTwo system was designed with your practice in mind. Excellent ROI potential — with low acquisition and disposable costs.

Positive outcomes – satisfied patients

- Treat acne scars & striae²
- Fractional skin resurfacing¹
- Treat facial wrinkles²

Sublime and Sublative Complementary Therapy

Wrinkles



Before Post 2 treatments
Photos courtesy of Ruthie Amir, M.D.
All photos are unretouched. Individual results may vary.

Wrinkles and Contouring



Before Post 1 treatment
Photos courtesy of Brian Biesman, M.D.
All photos are unretouched. Individual results may vary.

Sublime Wrinkles



Before Post 2 treatments
Photos courtesy of Ruthie Amir, M.D.
All photos are unretouched. Individual results may vary.

Sublative Striae



Before Post 3 treatments
Photos courtesy of Tsai, M.D.
All photos are unretouched. Individual results may vary.

eTwo combines its Sublime and Sublative applicators for aesthetic effect

The science of Sublime



Sublime uses infrared IR and bi-polar radio frequency (RF) energy to provide controlled deep dermal heating. This predictable heating of the deep dermis both stimulates new collagen growth and shortens aged collagen fibers, resulting in the immediate and long term improvement.

- Improvement of wrinkles and skin irregularities
- Smoother skin appearance
- Stimulates collagen production

The science of Sublative



The Sublative delivers fractionated bi-polar radio frequency (RF) energy via the matrix of electrodes of the disposable Sublative tip. This fractional delivery of RF energy causes micro-epidermal ablative injuries with wider spread dermal coagulation. This initiates a strong wound healing response in the dermis while also minimizing the damage to the epidermis. Less affected tissue surrounding the matrix points may be stimulated by heat which helps accelerate the healing process, leading to a smoother, rejuvenated appearance of the skin.

- Improvement in superficial skin tone and texture
- Treat facial wrinkles
- Use on ALL skin types

